

Transferring the Guardant DL[...]

07/16/2024 07:41:00

[FAQ Article Print](#)

| | | | |
|------------------|-------------------------------------|---------------------|---------------------|
| Category: | FAQ in English::Software Protection | Last update: | 06/26/2024 11:42:14 |
| State: | public (all) | | |

Problem (public)

Transfer the Guardant DL software key from one PC to another.

Solution (public)

The latest versions of the utility must be on the PC with the key and on the PC where you want to transfer the key.

The current version of the activation utility for Windows OS (x64 only):
[1]https://sd7.ascon.ru/Public/Utils/Guardant_SLK/kActivation/kActivation.zip

To transfer the software key: Create a transfer file on the PC where the key is located and apply this file to the PC where the key needs to be transferred.

1) Creating a transfer file .rehost Launch the activation application, select the TRANSFER KEY menu item. In the list of keys available for transfer, select the key to transfer and specify the path to save the transfer file .rehost.

2) Application of the transfer file on a PC with an Internet connection Launch the activation application, select the TRANSFER KEY menu item. Click Apply .rehost and specify the path to the .rehost transfer file.

[1] https://sd7.ascon.ru/Public/Utils/Guardant_SLK/kActivation/kActivation.zip